

OPTIONAL

Carbonic and Skin contact maceration

Available for the following models:

- | | | |
|-------------|----------------|--------------|
| - Velvet 23 | - Velvet 80 | - Velvet T43 |
| - Velvet 34 | - Velvet T11.5 | - Velvet T60 |
| - Velvet 40 | - Velvet T15 | - Velvet T75 |
| - Velvet 50 | - Velvet T26 | |
| - Velvet 65 | - Velvet T32 | |

Technical Description

This device allows to carry out the carbonic maceration and skin contact directly inside the pneumatic press. This device automates the homogenization of the grape, controls the pressure inside the tank and avoid the outdoor contact.

This system, thanks to a devoted software, controls the pressure inside the tank and rules gas relief that takes shape while the machine is stopped during the maceration phase. It is possible periodically to program the tank rotation in order to make the mass of grape uniform while macerating.